SOUND REACH H\F Shears

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
Rev. A.0


Reach Surgical, Inc. 120 Xinxing Road, West Zone, TEDA, 300462 Tianjin, China
EC REP MDSS GmbH Schiffgraben 41, 30175 Hannover, Germany

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en

Disposable Ultrasonic Scalpels Instructions

Please read all information carefully.

Failure to properly follow these instructions may lead to serious consequences.

Important: This package insert is designed to provide the instruction for use of **Disposable Ultrasonic Scalpels**.

Models include: CS1405P, CS2305P, CS3605P, CS4505P, CS1405H, CS2305H, CS3605H, CS4505H, CS1405F, CS2305F, CS3605F, and CS4505F.

It is not a reference for surgical techniques.

Chapter1 Product Overview

The operation manual is intended to help use the Disposable Ultrasonic Scalpels (hereinafter referred to as the "scalpel"), used as technical reference for non-surgical operations. The product is designed, inspected and manufactured for use in a single surgical operation. The instruments are EO sterile, single use device and single patient use Instrument used for dissection, grasping, coagulation, and cutting between the blade and clamp. If repeatedly used or reworked, the product may lose its functions, or may simultaneously cause patient's injury or cross contamination. Please never repeatedly use or rework the product.

Intended Use

This instrument is used to dissect soft-tissue for cutting and/or coagulating tissues.

Indications

This instrument is designed to be attached to an ultrasound surgical equipment and contact with a patient during vibration at high frequency in order to dissect soft-tissue cells for cutting and/or coagulating tissue in open surgeries or endoscopic surgeries. In general, pediatric, gynecologic, urologic, thoracic, and sealing and transection of lymphatic vessels.

Intended User

This instrument is intended to be used for healthcare professionals who use this instrument for surgical purposes.

Intended Use Environment

This instrument is intended to be used in a hospital.

Intended patient population:

Patients aged 3 and older who need surgery in which soft tissue incisions with bleeding control and minimal thermal injury are required.

Clinical benefit:

- Shorter operative time;
- Less intraoperative bleeding;
- Less thermal injury.

Contraindications

- * The instruments are not indicated for incising bone.
- * The instruments are not intended for contraceptive tubal occlusion

Patient target population

Patients who need surgery in which soft tissue incisions with bleeding control and minimal thermal injury are required.

Chapter2 Operation instructions

Scalpel consists of end-effector, rod and handle. The scalpel has multiple functions: its head can be used for spot homeostasis, its blunt side on one side can implement the best homeostasis function, simultaneously with tissue incision function, while the other sharp side functions rapid tissue incision and separation.

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Assembly

- The clamp is designed in its open position in the package; do not try to close it; do not apply too much force on the trigger.
- Align the screwed end of transducer with threaded hole in the handle along the direction shown in Fig. 2; hold transducer with one hand and hold the torque spanner of handle with the other hand, and turn clockwise till you can not turn.
- For type P**, close the clamp through the trigger; push the torque spanner along direction shown in Fig. 3 to the rotation knob of handle; turn the torque spanner clockwise till the set force is reached, and screw out the rod reversely.

For types H and F, the rotation knob has function equal to torque spanner; turn the rotation knob clockwise till the set force is reached, you will hear the click sound.

Caution: For type P, keep the clamp closed while the torque spanner is sleeved-in and the shaft is screwed-out.

Unloading

- For type P: Keep the clamp closed; sleeve the torque spanner on rod; hold the transducer with one hand, and screw off the transducer with the other hand anti-clockwise;** For types H and F: Turn the rotation knob anti-clockwise till out of the shaft.
- Hold the transducer with one hand while holding the rotation knob with the other hand, and screw the transducer out of the shaft.

Procedure for Use

Caution: Before the scalpel is inserted into the puncture outfit, the clamp must be closed, which can be implemented through the trigger.

- Insert the ultrasonic scalpel into a proper puncture outfit, or into a larger puncture outfit through an adapter;
 - After the ultrasonic scalpel is fully entered into coelom, open the clamp;
- Caution: Do not apply force on the trigger arbitrarily. It is not required to be inserted into the puncture outfit in the open surgeries.**
- The shaft can be turned by 360°. By rotation knob, adjust the scalpel end-effector to optimum position;
 - Place tissue to be cut into clamp and scalpel end-effector, and be sure that there are no other obstacles in the clamp mouth; length to be cut shall not exceed exposed scalpel end-effector length. Longer tissue to be cut may be incised several times;
 - Hold on trigger tightly, so as to close clamp jaws, in which tissue to be incised is contained.
 - After incision is completed, stop driving, open the clamp mouth, and remove the scalpel carefully. Check haemorrhage at stitch location. If there is any haemorrhage, manually suture and legate for haemostasis.
 - Close the clamp mouth and take the scalpel out of coelom. More system usage information can be referred to the operation manual of Ultrasound Surgical Equipment or Transducer.

Chapter 3 Warnings and Precautions

- The product can only be used for haemostasis of blood vessel of 3mm or smaller.
- If you fail to observe haemostasis situation at stapling location, do not use the product.
- The product is a single-use aseptis apparatus. Discard it after use. Do not sterilize and use it repeatedly.
- Before scalpel is inserted into or taken out of sleeve of puncture outfit, its clamp mouth shall be kept closed.
- During and after incision, haemostasis shall be checked; slight haemorrhage may, if any, be manually stitched for haemostasis.
- Tissue placed beyond exposed length of ultrasonic scalpel end-effector or too much force applied may cause apparatus failure and abnormal incision.
- Place tissue to be cut between clamp arm, and be sure that there are no other obstacles in the clamp mouth.
- After incision is finished, generator shall stop driving at the same time, and scalpel end-effector shall avoid contacting other tissues as possible.
- Blunt side is used for hemostasis at power level 3 or lower, or incision at power level 4 or higher. Sharp side of the scalpel is intended for incision at all power levels. Top of the scalpel is intended for spot hemostasis.
- Endoscopic surgery can only be performed by doctors with rich experience and familiar with endoscopic surgery. Before surgical operation, doctors shall consult related literatures so as to understand surgical technology, complication and risks.
- If apparatus and accessories from different manufacturers are used in a single operation, compatibility among them shall be verified, and insulation and grounding shall be checked for good condition.
- The product is sterilely packaged, only for use in a single operation. Discard it after use. Do not sterilize and use it repeatedly.
- Used products shall be placed in a specified recycle bin or refuse container. Do not litter to avoid environmental pollution.
- A notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to Reach Surgical, Inc. through Reachquality@reachsurgical.com. and the competent authority of the Member State in which the user and/or patient is established.
- Please use the product before its expiration date, and expired products are prohibited to be used.
- Scalpel of the Ultrasound Surgical Equipment must be stored in a room with ambient humidity not more than 80%, temperature -10°C ~ +55°C, atmospheric pressure 860hPa ~ 1060hPa, no acid or alkali or corrosive gas, good ventilation, no strong mechanical vibration or strong magnetic field or sunlight irradiation. After sterilization, the period of validity counted from sterilization is set to 5 years by following rules for storage.

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